



# PRODUCT CATALOG

SELECTED



Wuxi GALAK Chromatography Technology Co., Ltd




88th Meiliang Road  
Wuxi City, Jiangsu, 214092  
Mob: +86-1881-6200534  
Tel: +86-0510-85992929  
Fax: +86-0510-85104949  
[www.galaklc.com](http://www.galaklc.com)

GALAK Chromatography

**GALAK**  
CHROMATOGRAPHY

**Focus on Liquid Chromatography  
Since 2009**



We are continually developing new products to meet the needs of our customers.

For the complete GALAK product line including items not featured in this publication, please contact us directly or visit [www.galaklc.com](http://www.galaklc.com).

We are looking forward to helping you.

# GALAK Chromatography

Wuxi GALAK Chromatography Technology Co, Ltd (GALAK) is a technology-driven enterprise established in 2009 by experienced liquid chromatography experts who worked in the US and Japan. GALAK's headquarter is located in Wuxi Bio-Park with over 1,500 sq.m. of R&D center. Its factory is located in Zhuhai city, Guangdong Province with over 3,000 sq.m. equipped with state-of-art instruments and equipment. During 14 years, GALAK has established over 20 proprietary intellectual property rights on chromatography products and biochemical purification technologies.


GALAK products include HPLC prepacked columns, liquid chromatography packing materials, HPLC hardware, industrial purification systems for normal phase, reversed-phase, ion-exchange and affinity chromatography applications. For HPLC prepacked columns and liquid chromatography packing materials, Galaksil® silica-gel materials are silica products with C18, C8, CN, NH2, phenyl. Sepromax® materials are based on the PS-DVB matrix to isolate and purify peptides, proteins, polysaccharides, and antibodies. Vircap® materials are PS-DVB matrix to isolate and purify virus molecules.

With extended experiences in biochemical purification, GALAK helped over 50 companies and laboratories to design and improve their biochemical purification processes. The projects involved monoclonal antibodies purification, peptides isolation, protein isolation, nature chemical purification, enantiomers resolution and chromatography analysis.

"Innovation, Cooperation, Mutual benefits" are our philosophy.

GALAK is looking forward to work with you.

## Product Portfolio



### Prepacked Column

- Reversed-phase LC Column
- Normal-phase LC Column
- Absolut A Column



### Chromatography Resin

- Galaksil Silica-gel Resin
- Protein A Affinity Resin
- Sepromax Ion-exchange Resin
- VirCap® Perfusion Resin
- VirCap® Oligo dT(25) Affinity Resin
- VirCap® InertShell Core-Shell Resin



### Instrument & Accessories

- Column Packing System
- High-pressure Injection Pump
- Injection Loop
- Oligo Synthesis Column
- Glass Column
- HPLC Accessories

# OUR TECHNOLOGY

## Substrate Particle

Substrate particles build the foundation of the mechanical and chemical stability in packing materials.

## Sepromax® Polymer Particles

Sepromax® is a family of spherical divinylbenzene-styrenecopolymer (PS-DVB) particles designed for large-scale purification processes. With unique technologies, we precisely control their particle size, pore structure, pore size and surface area. Sepromax® particles have excellent mechanical properties and can withstand up to 10 MPa pressure. Their large pore sizes allow low mass transform of biomacromolecules.


Optional: average particle size: 20/50/70/150  $\mu\text{m}$ , pore size: 1000/2000/3000/5000A



## Viricap® Polymer Particles

Viricap® is also a family of spherical divinylbenzene-styrenecopolymer (PS-DVB) particles which designed for virus molecules purification processes. Viricap® particles are developed from Sepromax® particles, the difference is that the Viricap® particles have much larger pore structures and “through-pores”. These large through-pores allow part of mobile phase to flow through, quickly carrying biomolecules to smaller diffusive pores.

Optional: average particle size: 50/70/150  $\mu\text{m}$ , pore size: 3000/5000A



## Galaksil® Silica-gel Particles


Galaksil® is a well-established product series manufactured by innovative processes at industrial scales. It is a family of spherical, silica particles with tightly controlled particle size, pore structure and surface area.




## Ligand Bonding Technology

Functional groups or ligand on traditional packing materials are distributed along their surface. This limits their ability to effectively interact with bio-molecules. GALAK developed technologies to build “tentacle structure” on its ion-exchange and affinity packing materials. Flexible tentacle structures minimize the steric hindrance between functional groups or ligand and target molecules, improving the binding capability of the target molecules. Compared to traditional packing materials, GALAK’s PS-DVB particles show more effective capture and higher recovery.

**Innovation Ligand Tentacle**



**Traditional Ligand**



## Functional Groups

The functional group determines resin selectivity.

Common functional group classifications for silica-gel particles include reversed phase (RP), normal phase (NP), hydrophilic interaction chromatography (HILIC), ion-exchange (IEX), size exclusion chromatography (SEC), ion exclusion chromatography (ICE) and affinity chromatography (AC).

The functional group for PS-DVB particles include ion-exchange (IEX), and affinity chromatography (AC).

Affinity		
Recombinant Protein A	Sulfate Ester	dT-25mer
Cation-exchange		
Sulfonic Group S/SP	Strong Type	$-\text{SO}_3^{2-}$
Carboxymethyl CM	Weak Type	$-\text{COO}$
Anion-exchange		
Quaternary Ammonium Q	Strong Type	$-\text{N}^+(\text{CH}_3)_3$
Tertiary Amine D	Weak Type	$-\text{N}^+\text{H}(\text{CH}_3)_2$



## **We are GALAK Chromatography**

At GALAK, customer service isn't just a department—it's who we are. We are sales and quality assurance professionals, technical support specialists, engineers and chemists. We are here to support customers every step of the way. GALAK is committed to delivering innovative products that provide optimal performance.



## LC Prepacked Column

---

10	Reversed-phase LC Columns
18	Normal-phase LC Columns
21	Guard Column
22	Protein A Analysis Column

## Sorbents

---

22	<b>Silica-gel Packing Materials</b> Reversed-phase Chromatography Normal-phase Chromatography
29	<b>Affinity Chromatography</b> Protein A PS-DVB Affinity Resin Agarose Affinity Resin
37	<b>Ion-exchange Chromatography</b> PS-DVB IEX Resin Agarose IEX Resin
41	<b>Perfusion Chromatography</b> VirCap® AF Media VirCap® Oligo dT(25) Affinity Resin
49	<b>Multi-function Chromatography</b> VirCap® InertShell Core-Shell Resin

## Instruments & Parts

---

51	Packing System For HPLC Column
53	High-pressure Injection Pump
56	Glass Chromatography Column
61	Injection Loop
62	Oligo Synthesis Column
63	HPLC Column Hardware

# Prepacked Columns

Galaksil® prepacked columns are versatile HPLC columns based on the silica-gel for reversed-phase/normal phase chromatography. Galaksil® columns are made of spherical silica-gel particles which has low metal-ion content (<20 ppm in total), high specific surface area and high mechanical strength. With unique chemical bonding technique, our products have excellent stability and reproducibility. They can meet the highest requirements for analysis and preparative applications.

## Advantages

- Low silanol activity
- Uniform ligand binding
- Low metal content
- Narrow particle size




## Customized Columns







## Distribution of particle size for Galaksil® C18 5µm



## Repeated injection tests for Galaksil® C18 5µm



## National Institute of Standards and Technology (NIST) ERM 870 Test



### SRM 870 test

Column: Galaksil EF-C18H 150×4.6mmI.D

Eluent: 20mM Potassium Phosphate (pH7.0/ MeOH=20/80

Flow rate: 1.0mL/min

Det.:UV254nm

Temp.:23C

1.Uracil

2.Toluene

3.Ethylbenzene

4.Quinizarin

5.Amitriptyline

Galaksil® C18H can use in alkali environment with high pH CIP (Clean-in-Place) process. The isolation of toluene and ethylbenzene test shows the uniformities of binding ligands on the silica-gel substrate.

## Galaksil® Reversed-phase LC Columns

Galaksil® prepacked reversed-phase columns based on spherical silica-gel particles with very low metal-ion content (<20 ppm) in total, high specific surface area and high mechanical strength. Excellent base silica particles combined with unique chemical bonding technique ensures our products to have excellent stability and reproducibility. They can meet the highest requirements for analysis and preparative applications.

### Advantages:

- Low silanol activity
- Uniform ligand binding
- Low metal content
- Narrow particle size
- Excellent stability



Products	Particle Size	Pore Size	Surface Area	Carbon Content	pH Range
C18M	3/5/8/10 um	120Å	330m <sup>2</sup> /g	16%	2-8
C18H	5/8/10 um	120Å	330m <sup>2</sup> /g	20%	2-11
C18L	5/10 um	120Å	330m <sup>2</sup> /g	13%	2-8
C8	3/5/8/10 um	120Å	330m <sup>2</sup> /g	12%	2-8
C4 Bio	5/10um	300Å	100m <sup>2</sup> /g	3%	2-8
C8 Bio	5/10um	300Å	100m <sup>2</sup> /g	5%	2-8
C18 Bio	5/10 um	300Å	100m <sup>2</sup> /g	8%	2-8
Phenyl	3/5/10 um	140Å	300m <sup>2</sup> /g	8%	2-8



GALAK provide customized silica-gel sorbent, please contact us for details.

# Galaksil® C18M

## Parameters

Particle Size	Pore Size	Surface Area	Carbon Content	pH Range
3/5/8/10µm	120Å	330m <sup>2</sup> /g	16%	2-8

## Application



### Nucleotide

Column: C18M 5µm 4.6×150mm

Competitor ODS 5µm 4.6×150mm

### Mobile Phase:



phosphoric acid buffer / methyl alcohol

Flow Rate: 1ml/min

Wavelength: 254nm

Temp.: 25°C

1 5'-cytidylic acid; 2 5'-uridylic acid;  
3 5'-guanylic acid; 4 5'-inosinic acid;  
5 5'-adenylic acid



### Paraben

Column: EF-C18M 5µm

4.6×150mm

Competitor ODS 5µm

4.6×150mm


Mobile Phase: Water / methyl alcohol

Flow Rate: 1ml/min

Wavelength: 254nm

Temp.: 25°C

1 Methyl ester; 2 Ethyl ester;  
3 Propyl ester; 4 Butyl ester



### Water-soluble multivitamin

Column: C18M 5µm

4.6×150mm

### Mobile Phase:

phosphoric acid buffer / acetonitrile

Flow Rate: 1ml/min

Wavelength: 210nm

Temp.: 25°C



1 Pyridoxine;  
2 VB1;  
3 Nicotinamide;  
4 Folic acid;  
5 VB2

# Galaksil® C18H

## Parameters

Particle Size	Pore Size	Surface Area	Carbon Content	pH Range
5/8/10um	120Å	330m <sup>2</sup> /g	20%	2-11

## Application



### Ibuprofen/Benzene ketone



Column: EF-C18H 5µm 4.6×150mm  
Competitor 5µm 4.6×150mm

Mobile Phase:  
phosphoric acid buffer / acetonitrile

Flow Rate: 2ml/min

Wavelength: 214nm

Temp.: 30°C



### Polar/Nonpolar/ Neutral/Alkali Compounds

Column: EF-C18H 5µm 4.6×250mm  
Competitor 5µm 4.6×250mm

Mobile Phase:  
phosphoric acid buffer / methyl alcohol


Flow Rate: 1ml/min

Wavelength: 254nm

Temp.: 30°C

1 Uracil; 2 Butyl p-hydroxybenzoate;  
3 Propranolol; 4 Di-propyl ortho-phthalate;  
5 Naphthalene; 6 Acenaphthene;  
7 Amitriptyline

## The purification of EPA in fish oil




### EPA in fish oil

Column: C18H 8µm  
20×250mm

Sample: 90% EPA material

Finished sample

Purification: 99.7%






## Peptides Purification Test








Galaksil® UP-C18H and the word-leading competitive product in a peptides purification study. The results show that the Galaksil® UP-C18H is similar to the competitive product.

		Galaksil® C18	Competitor
Performance	Column Height (cm)	21.3	21.1
	Column Efficiency (TP)	70457	56935
Peptides	Injection Sample (g)	2.5	2.5
	Recovery (%)	89.3	90.0
	Purity(%)	95.5	95.5
	Freeze-dried product (g)	1.1302	1.1317

Galaksil® C18 8µm




Competitor



## Galaksi® C18L

### Parameters

Particle Size	Pore Size	Surface Area	Carbon Content	pH Range
5/10um	120Å	330m <sup>2</sup> /g	13%	2-8



#### Tripeptide (5ppm)

**Column:** C18L 5µm 4.6×250mm

**Mobile Phase:** 70/30 v/v Water/MeCN

**Injection:** 25µL

**Flow Rate:** 1ml/min

**Wavelength:** 220nm


**Temp.:** 25°C

## Galaksi® C8

### Parameters

Particle Size	Pore Size	Surface Area	Carbon Content	pH Range
3/5/8/10um	120Å	330m <sup>2</sup> /g	12%	2-8

### Application



#### Orlistat


**Column:** C8 5µm 4.6×250mm

**Mobile Phase:** water / EtOH

**Flow Rate:** 1ml/min

**Wavelength:** 203nm

**Temp.:** 25°C



#### Omeprazole enteric-coated tablets


**Column:** C8 5µm 4.6×250mm

**Mobile Phase:** water / EtOH

**Flow Rate:** 1ml/min

**Wavelength:** 203nm

**Temp.:** 25 °C



### Orlistat

**Column:** EP-C8 10 $\mu$ m 10 $\times$ 250mm

**Mobile Phase:** EtOH solution

**Flow Rate:** 4ml/min

**Wavelength:** 195nm

### Sample:

Dissolved raw material with methyl alcohol


**Concentration:** 50-60mg/ml

**Finished sample**

**Purification:** 99.8%

**Single impurity** < 0.1%

**Recovery:**  $\geq$ 90%



### Insulin

**Column:** C8 8 $\mu$ m 10 $\times$ 250mm

Time	A	B
0	85%	15%
5min	85%	15%
15min	64%	36%
225min	34%	66%

Galaxsil® C8	Cycle	Injection	Purification	P1	P1c	P2
	1	100ml	99.76%	0.21%	0.02%	0.01%
		50ml	99.74%	0.22%	0.02%	0.02%
	2	50ml	99.75%	0.22%	0.02%	0.01%
	3	50ml	99.74%	0.22%	0.02%	0.01%
	4	50ml	99.74%	0.22%	0.02%	0.01%
	5	50ml	99.76%	0.21%	0.02%	0.01%
	6	50ml	99.75%	0.22%	0.02%	0.02%
	7	50ml	99.76%	0.21%	0.02%	0.02%
	8	50ml	99.74%	0.22%	0.02%	0.01%
	9	50ml	99.74%	0.22%	0.02%	0.02%

## Galaksil® C4Bio

### Parameters

Particle Size	Pore Size	Surface Area	Carbon Content	pH Range
5/10µm	300Å	100m <sup>2</sup> /g	3%	2-8

## Galaksil® C8Bio


### Parameters

Particle Size	Pore Size	Surface Area	Carbon Content	pH Range
5/10µm	300Å	100m <sup>2</sup> /g	5%	2-8

## Galaksil® C18Bio

### Parameters

Particle Size	Pore Size	Surface Area	Carbon Content	pH Range
5/10µm	300Å	100m <sup>2</sup> /g	8%	2-8



#### Riboviron

**Column:** C18Bio, 5 µm 4.6×150 mm

**Mobile Phase:**

A) Na<sub>2</sub>SO<sub>4</sub>, pH2.5;

B) 40/60 v/v MeCN/Na<sub>2</sub>SO<sub>4</sub>, pH2.5

**Gradient:**

t (min)	%A	%B
0	100	0
15	100	0
25	87	13
35	87	13
50	0	100

**Flow Rate:** 1.0 mL/min

**Temperature:** 30°C

**Injection:** 10 µL

**Detection:** UV 220 nm

**Peaks:**


1. triazolinic acid;
2. Triazolamide;
3. Ribavirin acid;
4. Ribavirin;
5. Ribavirin 5 isomers;
6. Ribavirin methyl ester;
7. Ribavirin 5' - acetyl;
8. Ribavirin 5' - benzoyl



# Galaksil® Phenyl

## Parameters

Particle Size	Pore Size	Surface Area	Carbon Content	pH Range
5/10µm	140Å	300m <sup>2</sup> /g	8%	2-8



### Roflumilast

**Column:** Phenyl 5µm  
4.6×250mm

**Mobile Phase:** 60/40 v/v  
Water/MeCN

**Injection:** 10µL

**Flow Rate:** 1ml/min

**Wavelength:** 215nm

**Temp.:** 30°C

## Order Information

	2.1-50mm	2.1-150mm	4.6-50mm	4.6-150mm
EF-C18M 3µm	721-03012-002105	721-03012-002115	721-03012-004605	721-03012-004615

	4.6-150mm	4.6-250mm	10-250mm	20-250mm	30-250mm
EF-C18M 5µm	721-05012-004615	721-05012-004625	721-05012-010025	721-05012-020025	721-05012-030025
EF-C18H 5µm	722-05012-004615	722-05012-004625	722-05012-010025	722-05012-020025	722-05012-030025
EF-C18L 5µm	723-05012-004615	723-05012-004625	723-05012-010025	723-05012-020025	723-05012-030025
EF-C8 5µm	725-05012-004615	725-05012-004625	725-05012-010025	725-05012-020025	725-05012-030025
EF-C4Bio 5µm	730-05012-004615	730-05012-004625	730-05012-010025	730-05012-020025	730-05012-030025
EF-C8Bio 5µm	729-05012-004615	729-05012-004625	729-05012-010025	729-05012-020025	729-05012-030025
EF-C18Bio 5µm	728-05012-004615	728-05012-004625	728-05012-010025	728-05012-020025	728-05012-030025
EF-Phenyl 5µm	706-05012-004615	706-05012-004625	706-05012-010025	706-05012-020025	706-05012-030025

	4.6-250mm	10-250mm	20-250mm	30-250mm	50-250mm
EP-C18M 10µm	721-10012-004625	721-10012-010025	721-10012-020025	721-10012-030025	721-10012-050025
EP-C18H 10µm	722-10012-004625	722-10012-010025	722-10012-020025	722-10012-030025	722-10012-050025
EP-C8 10µm	725-10012-004625	725-10012-010025	725-10012-020025	725-10012-030025	725-10012-050025
EP-C4Bio 10µm	730-10012-004625	730-10012-010025	730-10012-020025	730-10012-030025	730-10012-050025
EP-C8Bio 10µm	729-10012-004625	729-10012-010025	729-10012-020025	729-10012-030025	729-10012-050025
EP-C18Bio 10µm	728-10012-004625	728-10012-010025	728-10012-020025	728-10012-030025	728-10012-050025

# Galaksil® Normal-phase LC Columns

Galaksil® prepacked normal –phase columns are made of spherical silica-gel particles with very low metal-ion content (<20 ppm in total) , high specific surface area and high mechanical strength. With chemical bonding technique, our products have excellent stability and reproducibility. They can meet the highest requirements for analysis and preparative applications.

## Advantages:

- Low silanol activity
- Uniform ligand binding
- Low metal content
- Narrow particle size
- Excellent stability


Products	Particle Size	Pore Size	Surface Area	Carbon Content	pH Range
SiO <sub>2</sub>	3/5/8/10 um	120Å	330m <sup>2</sup> /g	-	2-8
SiO <sub>2</sub> Bio	5 um	300Å	100m <sup>2</sup> /g	-	2-8
NH <sub>2</sub>	3/5/10 um	120Å	330m <sup>2</sup> /g	5%	2-8
CN	3/5/10 um	120Å	330m <sup>2</sup> /g	7%	2-8
Diol	5/10 um	120Å	330m <sup>2</sup> /g	8%	2-8
Amide	5/10 um	120Å	330m <sup>2</sup> /g	8%	2-8



## Galaksil® SiO<sub>2</sub>

### Parameters

Particle Size	Pore Size	Surface Area	pH Range
3/5/10µm	120Å	330m <sup>2</sup> /g	2-8



#### Maleic Maleic Fumaric Acid

**Column:** Galaksil SiO<sub>2</sub> 5µm  
4.6×250mm

**Mobile Phase:**  
N-hexane/THF/Trifluoroacetic  
acid = 650/350/1.2

**Injection:** 20µl

**Flow Rate:** 0.8ml/min

**Wavelength:** 255nm

**Temp.:** 30°C

## Galaksil® SiO<sub>2</sub> Bio



### Parameters

Particle Size	Pore Size	Surface Area	pH Range
3/5/10µm	120Å	330m <sup>2</sup> /g	2-8

## Galaksil® CN

### Parameters

Particle Size	Pore Size	Surface Area	Carbon Content	pH Range
3/5/10µm	120Å	330m <sup>2</sup> /g	7%	2-8



#### Benzalkonium Chloride

**Column:** Galaksil CN 5µm 4.6×150mm  
Competitor CN 5µm 4.6×150mm

**Mobile Phase:**  
phosphate buffer / acetonitrile

**Flow Rate:** 2.0ml/min


**Wavelength:** 214nm

**Temp.:** 35°C

## Galaksil® NH<sub>2</sub>

### Parameters

Particle Size	Pore Size	Surface Area	Carbon Content	pH Range
3/5/10µm	120Å	330m <sup>2</sup> /g	5%	2-8



#### Oligomaltose

**Column:** Galaksil NH<sub>2</sub> 5µm 4.6×150mm

**Mobile Phase:** water/ acetonitrile

**Flow Rate:** 1ml/min

**Detector:** RID

**Temp.:** 40°C

#### Peak

1 glucose; 2 maltose; 3 maltodextrin;

4 maltotetraose; 5 maltopentaose

## Galaksil® Diol

### Parameters

Particle Size	Pore Size	Surface Area	Carbon Content	pH Range
5/10µm	120Å	330m <sup>2</sup> /g	8%	2-8

## Order Information

	2.1-50mm	2.1-150mm	4.6-50mm	4.6-150mm
EF-SiO <sub>2</sub> 3µm	720-03012-002105	720-03012-002115	720-03012-004605	720-03012-004615
EF-NH <sub>2</sub> 3µm	705-03012-002105	705-03012-002115	705-03012-004605	705-03012-004615
EP-CN 3µm	704-03012-002105	704-03012-002115	704-03012-004605	704-03012-004615

	4.6-150mm	4.6-250mm	10-250mm	20-250mm	30-250mm
EF-SiO <sub>2</sub> 5µm	720-05012-004615	720-05012-004625	720-05012-010025	720-05012-020025	720-05012-030025
EF-NH <sub>2</sub> 5µm	705-05012-004615	705-05012-004625	705-05012-010025	705-05012-020025	705-05012-030025
EF-CN 5µm	704-05012-004615	704-05012-004625	704-05012-010025	704-05012-020025	704-05012-030025
EF-Phenyl 5µm	706-05012-004615	706-05012-004625	706-05012-010025	706-05012-020025	706-05012-030025
EF-Diol 5µm	707-05012-004615	707-05012-004625	707-05012-010025	707-05012-020025	707-05012-030025

	4.6-250mm	10-250mm	20-250mm	30-250mm	50-250mm
EP-SiO <sub>2</sub> 10µm	720-10012-004625	720-10012-010025	720-10012-020025	720-10012-030025	720-10012-050025




## Guard Columns

Cartridge + Holder

Size: 4.6-10mm, 10-10mm, 20-10mm

Packing material: matched with prepacked columns



## Precolumns

Size: 4.6-50mm, 10-30mm, 10-50mm, 20-30mm, 20-50mm, 30-50mm, 50-50mm

Packing material: matched with prepacked columns



# AbSolut® A Column

AbSolut® A column is designed for fast analysis of monoclonal antibody (mAb) concentration (titer) with protein A affinity chromatography. Alkali resistant recombinant Protein A (rProtein A) ligand used in this product has specific binding ability to the Fc region of immunoglobulins. The matrix of AbSolut® A is PS-DVB (Polystyrene Divinylbenzene) particles, which are highly cross-linked for enhanced mechanical stability and particle strength. Compared to agarose base, hydrophilic PS-DVB particles have higher pressure stability, dynamic binding capacity (DBC) and longer lifetime. Hence, AbSolut® A is an excellent choice for mAbs titer analysis.



## Advantages



- Direct use on HPLC instruments
- High dynamic binding capacity, quick mass transfer
- Minimum nonspecific absorption, accurate determination
- Fast analysis cycle time: 2–5 minutes
- Satisfactory linearity in wide concentration range: 0.02-10 mg/ml
- Long lifetime
- Alkali resistance: 0.1-0.5 M NaOH cleaning conditions

## Parameter

	AbSolut® A	AbSolut® A Plus
Column Size	2.1mm ID × 30mm L; 4.6mm ID × 50mm L	
Column Tube Material	316L Stainless steel, PEEK	
Support Matrix	Polystyrene Divinylbenzene (PS-DVB)	
Ligand	Recombinant Protein A	
Particle Size	30µm	20µm
Shipping Solution	0.02 M sodium phosphate, pH 7.0, 0.02% sodium azide	
pH range	pH 2-10	
Maximum Pressure	1000 psi	
Cleaning Agents	0.1-0.5M NaOH	
Cycle Time	2-5 minutes	
Temperature Stability	4-40 °C	

## Excellent Linearity


Quantitative analysis for antibody fermentation broth by Absolut® A.



**Column:** Absolut A 2×30mm  
**Eluent A:** 20mM PB, 150mM NaCl, pH7.4  
**Eluent B:** 0.1%HCl, 150mM NaCl  
**Gradient:** 0% B for 1.0 min, 100% B for 2.0 min, 0% B for 2.0 min  
**Flow rate:** 1ml/min  
**Sample:** mAb


## Flow Rate and Pressure

The operating flow rate is 0.5-3 ml/min as recommended for HPLC system.






**Column:** AbSolut A, 2.0×30mm  
**Eluent A:** 20mM PB, 150mM NaCl, pH7.4  
**Eluent B:** 0.1%HCl, 150mM NaCl  
**Temp:** 25 °C  
**System:** Waters 1525 pump

## Long Lifetime




<b>Column</b>	AbSolut A, 2×30 mm
<b>Eluent A</b>	50 mM Sodium Phosphate, 150 mM NaCl, pH 7.0
<b>Eluent B</b>	0.1% HCl, 150 mM NaCl, pH 1.9
<b>Flow Rate</b>	2.0 ml/min
<b>Gradient</b>	0% B for 0.2 min, 100% B for 0.60 min, 0% B for 1.20 min
<b>Temperature</b>	25°C
<b>Detection</b>	280 nm
<b>Injection volume</b>	10 uL
<b>Sample</b>	hIgG, 1 mg/mL

Statistical analysis of data demonstrates




## Flexible Choice of Flow Rate

The ratio of bounded and unbound IgG has almost no effect on the flow rate.



Normally, the flow rate is 1ml/min for 5 min analysis. Large samples, 2ml/min for 2 min analysis.




Flow rate ml/min	Total Area	Unbound Area	Unbound Relative Area %	IgG Area	IgG Relative Area %
0.5	1459568	145807	9.99	1313761	90.01
1.0	743661	75069	10.09	668592	89.91
1.5	492377	49715	10.01	442662	89.90
2.0	376354	39877	10.06	336477	89.40
2.5	322735	32984	10.22	289751	89.78



## Stability Test



Performance test for 10 different batches.

No.	RT (min)	Peak Area	Peak Height	TP	As
1	2.652	537586	190057	29507	1.10
2	2.641	536434	187236	26529	1.21
3	2.602	533688	186841	27349	1.12
4	2.599	531408	188244	29147	1.05
5	2.622	534911	187224	26901	0.98
6	2.647	540382	188746	26862	1.19
7	2.626	531906	188743	27855	1.08
8	2.628	540015	189618	28034	1.11
9	2.610	541372	188711	26567	1.16
10	2.623	527072	185477	26420	1.20





## Application Cases

Fc Fusion Protein sample



Monoclonal antibody sample



## Order Information




Particle Size	Column Size	Stainless Steel Column	PEEK Column
Absolut A	2.1* 30mm	1107-30-02003S	1107-30-02003P
	4.6*50mm	1107-30-04605S	1107-30-04605P
Absolut A PLUS	2.1* 30mm	1107-20-02003S	1107-20-02003P
	4.6*50mm	1107-20-04605S	1107-20-04605P

# Sorbents

With medium and large production lines, GALAK has the capacity for 300kg/liter to 800kg/liter per month of silica-gel packing materials and PS-DVB particles. We can meet production demand from pharmaceutical factories to laboratories. To satisfy customers needs, we also provide customized particles (particle size, pore size, ligand) based on specific technical requirements.

## Products:

- Galaksil silica-gel packing materials
- Sepromax® A50 for protein A affinity chromatography
- Sepromax® Polystyrene/divinylbenzene (PS-DVB) resin (Antibodies over 160kDa, other big protein molecules)
- VirCap® Polystyrene/divinylbenzene (PS-DVB) resin (capsule membrane virus and other virus for vaccines)
- VirCap® Oligo dT(25) Affinity Resin (mRNA purification)



# Reversed-phase Chromatography

## Analysis Silica-gel Packing Materials

Products	Particle Size	Pore Size	Surface Area	Carbon Content	pH Range	End-caped
EF-C18	3 um	120Å	330m <sup>2</sup> /g	16%	2-8	Yes
EF-C18M	5 um	120Å	330m <sup>2</sup> /g	16%	2-8	Yes
EF-C18M	8 um	120Å	330m <sup>2</sup> /g	16%	2-8	Yes
EF-C18M	10 um	120Å	330m <sup>2</sup> /g	16%	2-8	Yes
EF-C18H	5 um	120Å	330m <sup>2</sup> /g	20%	2-11	Yes
EF-C18H	10 um	120Å	330m <sup>2</sup> /g	20%	2-11	Yes
EF-C18L	5 um	120Å	330m <sup>2</sup> /g	13%	2-8	Yes
EF-C8	8 um	120Å	330m <sup>2</sup> /g	12%	2-8	Yes
EF-C8	10 um	120Å	330m <sup>2</sup> /g	12%	2-8	Yes
EF-C4 Bio	5 um	300Å	100m <sup>2</sup> /g	3%	2-8	Yes
EF-C4 Bio	10 um	300Å	100m <sup>2</sup> /g	3%	2-8	Yes
EF-C8 Bio	5 um	300Å	100m <sup>2</sup> /g	5%	2-8	Yes
EF-C8 Bio	10 um	300Å	100m <sup>2</sup> /g	5%	2-8	Yes
EF-C18 Bio	5 um	300Å	100m <sup>2</sup> /g	8%	2-8	Yes
EF-C18 Bio	10 um	300Å	100m <sup>2</sup> /g	8%	2-8	Yes
EF-Phenyl	5 um	120Å	330m <sup>2</sup> /g	8%	2-8	Yes
EF-Phenyl	10 um	120Å	330m <sup>2</sup> /g	8%	2-8	Yes

## Preparative Silica-gel Packing Materials

Products	Particle Size	Pore Size	Surface Area	Carbon Content	pH Range	End-caped
EP-C18	20 um	120Å	330m <sup>2</sup> /g	16%	2-8	Yes
EP-C18	30 um	100Å	330m <sup>2</sup> /g	16%	2-8	Yes
EP-C18	50 um	120Å	330m <sup>2</sup> /g	16%	2-8	Yes
EP-C18	75 um	120Å	330m <sup>2</sup> /g	16%	2-8	Yes
EP-C8	20 um	120Å	330m <sup>2</sup> /g	12%	2-8	Yes
EP-C8	30 um	100Å	330m <sup>2</sup> /g	12%	2-8	Yes
EP-C8	50 um	120Å	330m <sup>2</sup> /g	12%	2-8	Yes
EP-C8	75 um	120Å	330m <sup>2</sup> /g	12%	2-8	Yes

# Normal-phase Chromatography

## Analysis Silica-gel Packing Materials

Products	Particle Size	Pore Size	Surface Area	Carbon Content	pH Range	End-capped
EF-SiO <sub>2</sub>	3 µm	120Å	330m <sup>2</sup> /g	-	2-8	-
EF-SiO <sub>2</sub>	5 µm	120Å	330m <sup>2</sup> /g	-	2-8	-
EF-SiO <sub>2</sub>	7 µm	120Å	330m <sup>2</sup> /g	-	2-8	-
EF-SiO <sub>2</sub>	10 µm	120Å	330m <sup>2</sup> /g	-	2-8	-
EF-SiO <sub>2</sub>	5 µm	80Å	330m <sup>2</sup> /g	-	2-8	-
EF-SiO <sub>2</sub>	5 µm	100Å	330m <sup>2</sup> /g	-	2-8	-
EF-SiO <sub>2</sub>	10 µm	200Å	330m <sup>2</sup> /g	-	2-8	-
EF-SiO <sub>2</sub>	5 µm	300Å	330m <sup>2</sup> /g	-	2-8	-
EF-SiO <sub>2</sub>	10 µm	300Å	330m <sup>2</sup> /g	-	2-8	-
EF-NH <sub>2</sub>	5 µm	120Å	330m <sup>2</sup> /g	5%	2-8	Yes
EF-NH <sub>2</sub>	10 µm	120Å	330m <sup>2</sup> /g	5%	2-8	Yes
EF-CN	5 µm	120Å	330m <sup>2</sup> /g	7%	2-8	Yes
EF-CN	10 µm	120Å	330m <sup>2</sup> /g	7%	2-8	Yes
EF-Diol	5 µm	120Å	330m <sup>2</sup> /g	8%	2-8	Yes
EF-Diol	10 µm	120Å	330m <sup>2</sup> /g	8%	2-8	Yes


## Preparative Silica-gel Packing Materials

Products	Particle Size	Pore Size	Surface Area	Carbon Content	pH Range	End-capped
EF-SiO <sub>2</sub>	20 µm	120Å	330m <sup>2</sup> /g	-	2-8	-
EF-SiO <sub>2</sub>	30 µm	100Å	330m <sup>2</sup> /g	-	2-8	-
EF-SiO <sub>2</sub>	50 µm	120Å	330m <sup>2</sup> /g	-	2-8	-
EF-SiO <sub>2</sub>	75 µm	120Å	330m <sup>2</sup> /g	-	2-8	-

## Sepromax® A50

Sepromax® A50 is designed for analysis and purification of monoclonal antibodies (mAbs). Compared to traditional agarose media Sepromax® A50 has the advantages of high dynamic binding capacity (DBC), long life time, and less shedding of ligand. NaOH (0.1-0.5M) can be used for clean-in-Place (CIP).

The ligand of Sepromax® A50 is recombinant protein A (rProtein A) immobilized on the surface of macro-porous PS-DVB microsphere substrate. The rProtein A has better alkali-resisting ability that ensures stability in high pH conditions. With our hydrophilic treatment and coupling technology, we eliminated non-specific binding PS-DVB surface. Hence, Sepromax® A50 is extremely useful for purification process of monoclonal antibodies.



### Advantages :


- High rigidity, low backpressure, suitable for small-scale and large-scale mAb purification
- Outstanding high dynamic binding capacity at low residence time
- Excellent alkali resistance, 0.1-0.5 M NaOH for CIP
- Long lifetime, low ligand leakage
- High batch stability

### Parameter

<b>Support Matrix</b>	Poly(styrene/divinylbenzene) (PS-DVB)
<b>Ligand</b>	Recombinant Protein A
<b>Ave. Particle size</b>	50µm
<b>Dynamic Binding Capacity (DBC)</b>	Approx. 40 mg human IgG/ml media (Determined at 10% breakthrough by frontal analysis at a mobile phase velocity of 500 cm/h in a column with a bed height of 5 cm, Residence time 0.6 min)
<b>Shrinkage/Swelling</b>	< 1% from 1-100% organic solvent
<b>pH range (Long term)</b>	pH 2-10
<b>Maximum Operating Pressure</b>	1500 psi (100 bar / 10 MPa)
<b>Cleaning Agents</b>	0.1-0.5M NaOH
<b>Temperature Stability</b>	4-40 °C
<b>Delivery Conditions</b>	20% ethanol (2-8°C)

## DBC vs. Linear Flow Rate Curve


DBC of Sepromax<sup>®</sup> A50 do not decrease with the increased flow rate.



Condition	
System:	AKTA <sup>®</sup> Purifier10
Buffer:	20 mM PB, 0.15 M NaCl pH 7.4
Eluant:	0.1 M Gly-HCl, pH 2.5
CIP Solvent:	0.5 M NaOH
Sample:	1.0 mg/ml hlgG

## DBC vs. Residence Time Curve

Column efficiency of Sepromax<sup>®</sup> A50 is much higher compared to competitive products.



DBC=40 mg/ml, RT=1.3 min (Sepromax<sup>®</sup> A50),  
RT=2.4min (competitive product).

Column height=10cm, F=460cm/h (Sepromax<sup>®</sup> A50),  
F=250 cm/h (competitive product).

F=300 cm/h, H=6.5cm (Sepromax<sup>®</sup> A50),  
H=12 cm (competitive product).

$$RT = \frac{H \times 60}{F}$$

**Note:** Because of the difference of columns, residence time (RT) is more suitable for DBC description. F is for flow rate, H is for column height.

## DBC in High Sample Concentration

DBC of Sepromax<sup>®</sup> A50 is about 20-50% higher than competitive product with high IgG concentration injection under 2 minutes residence time (RT).


Test Condition	
Sample:	hlgG
Column:	7 mm I.D. x 2.5 cm (1mL)
Condition:	0.02 mol/L Na <sub>3</sub> PO <sub>4</sub> buffer (pH 7.4) + 0.15 mol/L NaCl
DBC:	Base on breakthrough curve (allow 5% leakage)

Residence time (min)	Flow rate (ml/min)	DBC @ 5% BT			
		IgG conc. at 5g/L		IgG conc. at 10g/L	
		Sepromax <sup>®</sup> A50	Competitor	Sepromax <sup>®</sup> A50	Competitor
1.92	0.5	30.3	27.4	30.3	24.5
0.96	1.0	26.8	15.3	26.7	14.5
0.64	1.5	23.2	9.8	23.4	10.8




## Pressure vs. Flow Rate Curve

The back pressure of Sepromax® A50 is less than competitive product. Sepromax® A50 is more suitable for industrial purification process.




## Protein A Ligand Leakage Test (ELISA)

Samples: 20mg hlgG/ml-resin, 9mg mAb/ml-resin; ELISA test: Cygnus F400 kit




Protein A concentration standard curve



## NaOH Clean-in-Place Test




After 120 cycles of 0.5M NaOH CIP, the relative DBC is still stay in high level.




## Alkali Resistance Test (150 Cycle Lifetime)

Do CIP after purification process, measure DBC each 18 injection cycles.

Step	solution	CV
Equilibration	20 mM PB, 0.15M NaCl, pH7.4	5CV
Loading	1mg/ml hlgG+ BSA/Lysozyme	RT=0.6min, 50%DBC
Washing	20 mM PB, 0.15M NaCl, pH7.4	5CV
Elution	0.1M Gly-HCl, pH3.0	5CV
CIP	0.1M NaOH	4CV, 15min






- Lane**
- 1: Marker
  - 2: Sample
  - 3: Flow through
  - 4: Run 4
  - 5: Run 21
  - 6: Run 39
  - 7: Run 57
  - 8: Run 75
  - 9: Marker
  - 10: Run 95
  - 11: Run 111

### Load-Eluent Recovery

With the rProtein A ligand fully functional, Sepromax<sup>®</sup> A50 delivers high recovery of purified antibodies.



Sample: 1.0mg/ml  $\gamma$ -globulin

	Load ( $\mu$ g)	Eluent ( $\mu$ g)	Recovery (%)
Sepromax <sup>®</sup> A50	37.65	33.16	88.07
Sepromax <sup>®</sup> A40 Plus	36.01	31.97	88.61
Competitor	33.96	30.09	88.81




### mAb Purification Test

Sample: monoclonal antibody from murine



- Lane**
- 1 Ascites treatment fluid (reduction)
  - 2 break through (reduction)
  - 3 break through (reduction)
  - 4 eluent (reduction)
  - 5 marker

### Fc protein Purification Test



<b>Media</b>	Sepromax A 50 affinity media
<b>Column</b>	1.0 × 2.5cm, column volume 2mL
<b>Sample</b>	Fc protein Xsupernatant(2.0g/L)
<b>Loading buffer</b>	10mM PB+0.2M NaCl, pH 7.5
<b>Elution buffer</b>	20mM Sodium Citrate+0.2M NaCl, pH3.7
<b>Flow rate</b>	115 cm/h


## Impurity Removal Test (HCP&DNA)

In the production of mAb's for pharmaceutical applications, residue of host protein (HCP) and DNA are an important indicators of quality. Protein A affinity chromatography is an efficient method to remove these residual impurities.


<b>Media</b>	1.0 × 2.5cm , column volume 2mL Sepromax A50
<b>Sample</b>	Fc protein Xsupernatant(2.0g/L)
<b>Loading buffer</b>	10mM PB+0.2M NaCl, pH 7.5
<b>Elution buffer</b>	20mM Sodium Citrate+0.2M NaCl, pH3.7
<b>Flow rate</b>	115 cm/h

### Result:

HCP (ng/mg)	Sepromax® A50	Competitor
supernatant	1754.3	1716.5
rPA Eluate	1.9	4.5
Reduction	$9.2 \times 10^2$	$3.8 \times 10^2$



Standard Curve of HCP



Standard HCP µg


Sample	Concentration	Wells	Values	Mean Value	Std. Dev.	CV%
HC01	0.000	A4	0.052	0.056	0.004	8.0
		A6	0.059			
HC02	0.700	B4	0.112	0.112	0.000	0.3
		B5	0.113			
HC03	2.000	C4	0.206	0.212	0.008	4.0
		C5	0.218			
HC04	8.000	D4	0.588	0.589	0.001	0.3
		D5	0.590			
HC05	25.000	E4	1.190	1.192	0.003	0.2
		E5	1.194			
HC06	75.000	F4	1.751	1.734	0.023	1.3
		F5	1.718			

Smallest standard value: 0.056  
Largest standard value: 1.734

## DNA Removal Test

DNA test kit: AB (4413713)

QPCR - using magnetic beads extracted DNA from sample. Prepare PCR reaction mixture with DNA extraction solution and standard solution. Using Bio-rad real-time PCR for reaction and fluorescence assay.



rDNA (pg/mg)	Sepromax® A40 plus G	Competitor A
Cell supernatant	$1.5 \times 10^7$	$1.6 \times 10^7$
rProtein A eluate	154.4	431.3
Reduction	$9.7 \times 10^4$	$3.7 \times 10^4$

## Regulatory Support Files, RSF

All regulatory support documents based on FDA reporting requirements that can assist customers in process development, validation and preparation of SOPs.

# GLK-gel Agarose Media

GLK-gel Agarose media offer high specificity and selectivity for biomolecular separations and purifications. Affinity separation can often remove contaminants difficult to eliminate using other chromatographic procedures. Purifications up to several orders of magnitude can be achieved in a single step.

## Advantages

- Stable bonding
- Low ligand leaching
- NaOH CIP

	Pr A 4FF	Pr G 4FF	Pr SupA
<b>Substrate</b>	Cross-linked agarose		
<b>Ligand</b>	rProtein A	rProtein G	Super alkaline resistant Protein A
<b>Particle Size</b>	90µm (45-165µm)		~60 µm
<b>Capacity (DBC)</b>	20mg hlgG/ml	25mg hlgG/ml	> 80 mg hlgG/ml
<b>pH Stability</b>	2-10 (Short) 3-9 (Long)		3-10
<b>Max. Pressure</b>	0.3MPa		
<b>Flow Rate</b>	300cm/h	300cm/h	100-500 cm/h
<b>Storage</b>	4-8 °C, 20% EtOH		1XPBS with 20% ethanol, 2-8°C

## Purification of IgG in human serum

**Sample:** 5ml human serum with five times dilution (different buffers)


**Column:** HT01 1.0ml Protein G 4FF

**Balance:** A 0.02 M PB pH7.0;



B 0.02M PB, 0. 3M NaCl pH 7.0

**Elution:** 0.1 M Glycine-HCl pH2.7

**Flow Rate:** 0.25m/min (sampling), 1ml/min



## Protein Purification



# GLK-gel Ni Affinity Media

GLK-gel Ni affinity media are a nickel metal chelating chromatography media with IDA/NTA/TED ion high cross-linked agarose. GLK gel Ni Affinity Media have advantages of excellent stability, biocompatibility, solvent compatibility, large capacity, good selectivity, high resolution natural generation and low cost.

## GLK-gel Ni IDA

<b>Substrate</b>	6% high cross-linked agarose
<b>Particle Size</b>	90µm (45-165µm)
<b>Binding Capacity</b>	Approx. 45 mg His (tag protein)/ml media
<b>pH Stability</b>	3-12 (Working); 2-14 (Cleaning)
<b>Max. Pressure</b>	0.3MPa
<b>Chemical Stability</b>	Common aqueous solutions and buffers. Avoid chelating agents (EDTA, EGTA) and reducing agents (DTT, DTE)
<b>Storage</b>	4-15 °C, 20°C EtOH

## Application Case

Column: 1ml

Sample: E. coli cracking supernatant (His tag protein)

Equilibrium liquid: 0.02MPB、0.5M NaCl, pH 7.4


Elution: 0.02MPB, 0. 5M NaCl, Imidazole, pH 7.4

Flow Rate: 1ml/min

Sample:

1. Original;
2. Breakthrough;
3. Elution(4%B);
4. Elution(100%B);
5. Elution(100%B);
7. Original;
8. Breakthrough;
9. Elution(4%B);
10. Elution(100%B)

No imidazole in 1-5. 0.02M imidazole in 7-10.



## GLK-gel Ni TED

Tolerance of higher reducing agents and chelating agents, eukaryotic secreted expression of His tag protein can loading without prior treatment, maximum protect the activity of protein.

Direct use NaOH for cleaning without nickel removal, reduce cleaning time.

Lower nickel shedding, no need for repeated regeneration.

<b>Substrate</b>	6% high cross-linked agarose
<b>Particle Size</b>	90µm (45-165µm)
<b>Binding Capacity</b>	Approx. 20 mg His (tag protein)/ml media
<b>pH Stability</b>	3-12 (Working); 2-14 (Cleaning)
<b>Max. Pressure</b>	0.3MPa
<b>Chemical Stability</b>	0.01M hydrochloric acid; 0.01M sodium hydroxide (one week); 20mM EDTA; 10mM DTT; 1M sodium hydroxide; 8M urea; 100mM EDTA; 0.5m imidazole (2 hours); 6M guanidine hydrochloride (24 hours); 30% isopropanol (20 min)
<b>Storage</b>	4-15 °C, 20°C EtOH



# Sepromax<sup>®</sup> IEX Media

## Alkali-resistance Type

Sepromax<sup>®</sup> ion-exchange media are based on large pore polystyrene/divinylbenzene (PS-DVB) particles. It has excellent mechanical property and can withstand pressures up to 10 MPa. Their 1000Å pore size allows low mass transform of biomacromolecules. These particles have been modified with GALAK unique coating technology and become completely hydrophilic.




### Advantages:

- Rigid particles, low backpressure, suitable for large-scale purification processes.
- High flow rate, high loading capacity, high purification efficiency.
- Excellent chemical stability, alkali stable under CIP and long lifetime.

## Sepromax<sup>®</sup> Polymer Particles



Sepromax<sup>®</sup> is a family of spherical divinylbenzene-styrenecopolymer (PS-DVB) particles designed for large-scale purification processes. With unique technologies, we precisely control their particle size, pore structure, pore size and surface area. Sepromax<sup>®</sup> particles have excellent mechanical properties and can withstand up to 10 MPa pressure. Their large pore sizes allow low mass transform of biomacromolecules.

Optional: average particle size: 20/50/70/150  $\mu\text{m}$ , pore size: 1000/2000/3000/5000Å



## Sepromax<sup>®</sup> “Tentacle”

Sepromax<sup>®</sup> IEX media carry “tentacle” surface structures. Functional groups are covalently bonded on the surface in the form of linear polymer chains. This structure enables macromolecules such as antibodies, viruses and plasmids to interact more effectively to the functional groups of the media, increasing the binding capacity significantly. “Tentacle” structure also effectively reduces the non-specific interaction between biomolecules and media, thus improving the recovery of target molecules.



## Parameter

	Sepromax <sup>®</sup> S50	Sepromax <sup>®</sup> CM50	Sepromax <sup>®</sup> Q50	Sepromax <sup>®</sup> D50
<b>Substrate</b>	Rigid, PS-DVB microspheres			
<b>Particle Size</b>	50um (35-75µm)			
<b>Ligand</b>	-SO <sup>3-</sup>	-COO-	-N <sup>+</sup> (CH <sub>3</sub> ) <sub>3</sub>	-N <sup>+</sup> H(CH <sub>3</sub> ) <sub>2</sub>
<b>pH Range</b>	2-12	6-12	2-12	2-9
<b>Dynamic Capacity</b>	60mg Lysozyme/ml	80mg Lysozyme/ml	100mg BSA/ml	100mg BSA/ml
<b>Max Pressure</b>	1500 psi (100 bar or 10 MPa)			
<b>pH Stability</b>	1-14			
<b>Storage</b>	20% EtOH, 4-30°C			


\* DBC (Dynamic Binding Capacity): frontal analysis @ 10%, 300cm/h, 5cm column height

## Customization Product

	Particle Size	Ligand	Dynamic Capacity	Application
Sepromax <sup>®</sup> SS50	50um (35-75µm)	-SO <sup>3-</sup>	80 mg Lysozyme/ml	Small molecular protein, <100,000kDa
Sepromax <sup>®</sup> S50			60mg Lysozyme/ml	Big molecular protein, >100,000kDa

## High Loading Capacity


Sepromax<sup>®</sup> S50 has excellent binding capacity under high linear velocity. This leads smaller column size and faster cycle time.



## Break-through Curve


The capture efficiency of Sepromax<sup>®</sup> S50 was measured by the fronted break-through curves at 5% and 10%.

The break-through point of protein penetration curve of polysaccharide type media is relatively earlier.






## Excellent Stability

The dynamic binding capacity (DBC) of Sepromax<sup>®</sup> IEX will not decrease significantly with high linear flow rate. When flow rate increased from 100cm/h to 500cm/h, Sepromax<sup>®</sup> S50 can maintain 65% of its DBC at 100cm/h.





## Stability under CIP Condition

Clean-in-place (CIP) is a very important process in protein purification in biopharmaceutical industry. Sepromax® IEX media shows excellent chemical stability under harsh CIP conditions. In the experiment, 1M NaOH solution was selected to soak the four ion exchange media, and the loading was evaluated at regular intervals. After soaking for one year, the loading capacity of the four media did not decrease significantly.



## Monoclonal Antibody Purification



Purify of finished product: 99.5%

Column: 2.5cm Height, 5ml CV

Medium: Sepromax® S50

Sample: Mab, Protein A elution pool 150 ml

Buffer A: 20 mM NaAc-HAC, pH 6.0

Buffer B: 20 mM NaAc-HAC, pH 6.0+ 1 M NaCl


Flow rate: 2 ml/min

Gradient: 0%-100% B (20 CV)


System: AKTA explorer100

# VirCap® Perfusion Media

Perfusion chromatography media are porous polymer microspheres made by polymerizing styrene and divinylbenzene. The particles contain two types of pores: large through-pores, also known as convective pores or perfusion pores, with diameters between 600-800nm that run throughout the entire particle, and smaller interconnecting pores, or diffusion pores, with diameters between 50-150nm and depths usually less than 1µm, that link the large pores. The media have a porosity of approximately 50-60%, can withstand pressures up to 5MPa, and can be used at higher flow rates than traditional media such as agarose or dextran (with pressure tolerance lower than 0.5 MPa), resulting in greatly improved production efficiency.



Perfusion chromatography media is a rigid polymer microsphere coated with a unique hydrophilic polymer, covalently linked to various functional groups (ion exchange, affinity, hydrophobic, etc.), with a highly stable structure that is very suitable for the purification of biopharmaceuticals. The large particle size provides a good balance between resolution and operational back pressure. The dynamic load and resolution are less affected by the increase in linear velocity, making the separation step faster than traditional gel chromatography media.



## Characteristic

- **Large Pore Size**

1000-3000Å pore size, enable the diffusion and mass transfer for large biomolecules.

- **Particle Size**

35-85 micron particles, satisfy your purification processing requests.

- **Rigid Microspheres**

Maximum pressure is over 870psi(60 bar), excellent mechanical property.

- **Flexible Tentacles**

Higher recovery rate and target purity, with excellent combination and capture capability .

- **Harsh Clean-in-Place Condition (CIP)**

0.5-1M NaOH, organic solvent, high salt solvent.





- **Robust Chemical Stability**

VirCap® particles are rigid polymeric particles that are coated with a proprietary hydrophilic polymer onto which the various functional groups (ion exchange, affinity, etc.) are covalently attached.





## Rigid Microsphere With Large Pore Size

VirCap<sup>®</sup> particles large “through-pores”. These large through-pores allow part of mobile phase to flow through, quickly carrying biomolecules to smaller diffusive pores. The large through-pores reduce diffusion rate of biomolecules and enhance interaction between biomolecules and functional groups on the surface. Consequently, mass transfer barriers are lowered, and flow rate can be increased without compromising capacity or resolution.




## High Pressure Resistance




## Flexible Tentacles

Flexible tentacle structures minimize steric hindrance between functional groups and target molecules. It also improves the binding capability of the target material. Compared to traditional media, VirCap<sup>®</sup> media show more effective capture and higher recovery.

Tentacle ion chromatography media



Traditional ion chromatography media







## High Flow Rate, Low Backpressure

VirCap® media offer an excellent balance of resolution and operating backpressure.

Under recommended condition of mobile phase, VirCap® media exhibit almost no shrinking or swelling. The combination of through-pores and flexible tentacles ensure rapid diffusion of solute. It also reduces the barrier of mass transfer, and realizes high dynamic binding capacity (DBC) under the high flow rate.



Peak width at different flow rates (VirCap® media 3000Å)




Dynamic capability at different flow rates (VirCap® media 3000Å)

## Robust Chemical Stability

VirCap® media are highly cross-linked polymeric particles coated with a proprietary hydrophilic layer on which various functional groups (ion exchange, affinity, etc.) are covalently attached. The result is chemically stable product that is ideally suitable for large-scale biopharmaceutical separation.

Lot	RT (min)	Area	Height	TP	As
1	2.652	537586	190057	29507	1.10
2	2.641	536434	187236	26529	1.21
3	2.602	533688	186841	27349	1.12
4	2.599	531408	188244	29147	1.05
5	2.622	534911	187224	26901	0.98
6	2.647	540382	188746	26862	1.19
7	2.626	531906	188743	27855	1.08
8	2.628	540015	189618	28034	1.11
9	2.610	541372	188711	26567	1.16
10	2.623	527072	185477	26420	1.20



## VirCap® AF / Q Media

Virus purification often used in producing virus type vaccines, and also provides an important tool for the study of virus fine morphological structure. Isolation and purification of virus antigen protein are detailed studies of virus chemical composition and genetic material.

VirCap® AF is an affinity chromatography media designed for the capture and moderate purification stages of capsular virus purification. Specific adsorption of VirCap® AF media and target occurs by simulating the affinity between ligands and virus particles with capsular membranes. With unique high loading capacity, high flow rate and low back pressure, VirCap® AF reduces the process cycle time and increases the yield, fully meeting the requirements of large-scale vaccine production processes. VirCap® Q on the other hand is a strong anionic exchange packing material that is capable of capturing virus type vaccine.

	VirCap® AF	VirCap® Q
<b>Substrate</b>	Hydrophilic PS-DVB (Polystyrene/divinylbenzene) Microspheres	
<b>Particle Size</b>	50um, 70um	
<b>Function Group</b>	Sulfate Ester	-CH <sub>2</sub> -CH <sub>2</sub> -CH <sub>2</sub> -N <sup>+</sup> (CH <sub>3</sub> ) <sub>3</sub>
<b>Dynamic Binding Capability</b>	lysozyme 30mg/ml	BSA >90 mg/ml
<b>Flow Rate</b>	1000cm/h (20°C, buffer solution viscosity same as water, pressure < 3 bar / 43.5psi, column bed height 20cm)	
<b>Column Bed Height</b>	20-40cm	
<b>pH Stability</b>	1-14	
<b>Working Temperature</b>	4-30°C	
<b>CIP Condition</b>	0.5-1M NaOH	
<b>Storage</b>	2-8°C 20% EtOH	

### VirCap® AF media Application

Viruses		Viral/Microbial Antigens
Rabies	Feline Calicivirus	Herpes Simplex gA and gB Glycoprotein Subunits
Influenza	Respiratory Syncytial Virus	Hepatitis B Surface Antigen
Japanese Encephalitis	Human Herpes Simplex	Filamentous Hemagglutinin from B. pertussis
Feline Leukemia	Human Measles	Leucocytosis Promoting Factor Hemagglutinin
Feline Herpes	Human Parainfluenza	


## One-step Porcine Pseudorabies Virus Purification

Porcine pseudorabies virus (PRV) causes fever, itchiness (except in pigs) and encephalomyelitis as the main symptoms in a variety of domestic and wild animals. Immunization is the main strategy for the prevention of pseudorabies, and a weakened vaccine with the Bartha-K61 strain is currently used in China.

The use of such inactivated virus vaccines is considered an effective way to prevent pseudorabies in pig farms, improving reproduction rate of sows and control piglet mortality.




The loading volume of VirCap® AF70 affinity chromatography is large (up to 5-10 column volumes). And it does not require concentration, which also avoids the loss of antigen from concentration and improves production efficiency. Therefore it is suitable to process scale-up.

The results of serum antibody detection showed that the antibody level after vaccination of purified vaccine - high dose group and medium dose group was close to that of commercial vaccine group, and the immunization effect was satisfactory.




### Advantages

1. Samples are pre-treated and directly sampled after VirCap one-step chromatography target yield greater than 70%; very low back pressure at higher flow rates.
2. Samples are loaded under neutral conditions with the vast majority of proteins, nucleases, HCP, endotoxins, DNA flow-through in the sample, and data provided by users indicate that the removal of miscellaneous proteins, HCP, nucleases, DNA, etc. is greater than 90%.
3. Mild adsorption and elution conditions, reducing downstream purification steps, sample pre-treatment only requires simple solid-liquid separation to remove most of the large solid particles.




1. Unprocessed original sample
2. Original sample centrifuged at 3000rpm
3. FT6
4. W6
5. ET6-1
6. FT7
7. W7
8. FT7-1
9. ET6-1 diluted 7 times
10. ET7-1 diluted 7 times
11. ET6-2 diluted 7 times
12. ET7-2 diluted 7 times
13. CIP6

## Gel Filtration Chromatography vs. Perfusion Chromatography



## SEC-HPLC Analysis



## Pilot Test


**Column:** VirCap AF 14.13L-CV, ID300mm L20cm.  
**Dilution buffer:** 20 mM MES pH 4-4.5  
**Loading buffer:** 20 mM MES pH 6.0-6.4  
**Post-sample Equilibration Buffer:** 20 mM MES pH 6.0-6.4  
**Elution buffer:** 50 mM Tris +0.5 M NaCl (pH~10)  
**Washing buffer:** 1M NaOH  
**Purification Flow Rate:** 140 L/h (200cm/hr)  
**Detection:** 280nm & 260nm

**Sample Pretreatment:** clarified liquid after centrifugation is filtered by 0.45um membrane, and 1X volume of dilution buffer is added to control sample conductance and pH (conductance range 7-10 mS/cm, pH range 5.9-6.3).

**Sample volume:** The recommended sample volume is 120-140 L (i.e. 8-10 CV).


**Elution:** After the sample loading is completed, first equilibrate the buffer to UV280 baseline and start elution with elution buffer to collect the elution peaks, and it is recommended that the collection termination point drops to 10% of the highest peak of the elution peak.

**Washing & Regeneration:** After elution is completed, wash 1-2CV with washing buffer, rest time 30-60min, equilibrate to neutral with pure water, and finally top sample/equilibration buffer 2-5CV, to be used.



# VirCap® Oligo dT(25) Affinity Resin

VirCap® Oligo dT(25) Affinity Resin is based on rigid, 50µm polymeric resin designed to isolate messenger RNA (mRNA). The resin backbone consists of crosslinked PS-DVB (polystyrene divinylbenzene).



The polyhydroxy surface coating provides low non-specific binding. The surface is functionalized with a linker and poly dT(25) functional group allowing capture of mRNA through H-bonding pairing with the mRNA polyA tail.

VirCap® Oligo dT(25) Affinity Resin provides efficient capture and easy release under standard mRNA purification conditions. It thereby decreases process development time and enhances productivity. In addition, the selective nature of this resin allows a reduction in plasmid DNA and other transcription mix components. The resin is also stable at elevated temperatures for the breakdown of undesired higher-order structures if required.

## Features

- Easy mRNA purification to separate non-poly A tail contaminants
- Simplified workflow helps to maximize efficiency, thereby reducing complexity of subsequent polish steps
- Excellent scalability
- Non-animal derived


*VirCap® Oligo ligands are manufactured using a synthetic manufacturing process that are free of animal components.*

## Specification

Characteristic	Description
<b>Support matrix</b>	cross-linked poly(styrene-divinylbenzene)
<b>Average particle size</b>	50 µm
<b>Average pore size</b>	200 nm
<b>Surface functionality</b>	poly(dT) 25mer with proprietary linker
<b>Ligand density</b>	0.3 µmol/ml
<b>Mechanical resistance</b>	70 bar (1,000 psi; 7 MPa)
<b>Thermal stability</b>	allows sample denaturing at 65°C if needed
<b>pH range</b>	2-13
<b>Ionic strength range</b>	0 to 5 M, all common salts
<b>Chemical resistance</b>	Common agents for mRNA purification, including 0.5 M NaOH, 2 M MgCl <sub>2</sub> , 20 mM EDTA. Water, 0 to 100% alcohol, acetonitrile, 2 M acetic acid, 1 M HCl, other common organic solvents
<b>Shipping solvent</b>	18-20% ethanol

## Capacity & Recovery



<b>Column size</b>	1ml column
<b>System</b>	AKTA purifier10
<b>Binding buffer</b>	10 mM Tris-HCl, 0.5 M NaCl, 1 mM EDTA, pH 7.4
<b>Wash buffer</b>	10 mM Tris-HCl, 300 mM NaCl, 1 mM EDTA, pH 7.4
<b>Elution buffer</b>	10 mM Tris-HCl, 1 mM EDTA, pH 7.4
<b>Regeneration</b>	Water
<b>CIP</b>	0.1M NaOH
<b>Sample</b>	mRNA
<b>Flow rate</b>	1ml/min



## Alkali Resistance Test


Method: After soaking the resin in 0.5 M sodium hydroxide for 48 hours, measure the change in binding capacity for poly A 40mer and mRNA.

Results: (VirCap® Oligo dT25 and 2 commercial Oligo dT resin)



## Batch Stability


Method: Five different batches of VirCap® Oligo dT(25) affinity resin were packed into columns and subjected to mRNA sample injection and elution program on HPLC. Chromatograms and peak areas of elution were compared.





# VirCap® InertShell Core-Shell Resin

VirCap® InertShell is designed with the core-shell technology. It is for purification of viruses and other large biomolecules. The core-shell technology allows for combining size exclusion separation with IEX chromatography. Viruses and other large biomolecules that are too large to penetrate the inert shell of the chromatography resins are collected in flow through fraction (FT mode). Contaminants (< Mr 700 000) on the other hand pass through the inert outer shell and bind to the ligands in the inner core.



VirCap® InertShell is made of polymethacrylate microspheres with octylamine ligand inside the pore as shown in Pic 1. The shell of the microspheres is neutral and hydrophilic, which has no biomolecule adsorption. The pore size of the shell (50-100 nm) is smaller than that of the core (200-500 nm). And the thickness of the shell is about 0.5-1.0µm. The shell prevents proteins (molecular weights greater than 700 kDa) from entering the core. In the chromatography process, large-size viruses or other large biomolecules cannot enter the microsphere core and they will breakthrough and be collected. The octylamine ligand in the core carries its dual functions of anion exchange and hydrophobicity, capturing protein molecules with molecular weight less than 700 kDa. VirCap® InertShell can effectively remove host cell proteins (HCP's), DNA fragment, endotoxin, albumin and other impurities.

## Specification

	VirCap® Inert Shell	Capto Core 700
<b>Matrix</b>	Polyacrylate	Highly cross-linked agarose
<b>Ligand</b>	Octylamine	Octylamine
<b>Average particle size</b>	50-150 µm	50-150 µm
<b>Density of ligand</b>	0.10-0.20 mmol/mL	0.04-0.085 mmol/mL
<b>Binding capacity<sup>1</sup></b>	6-12 mg BSA/mL resin	12 mg BSA/mL resin
<b>Operational pressure</b>	≤ 1.0 MPa	≤ 0.3 MPa
<b>Operational flow rate</b>	100-600 cm/h	100-600 cm/h
<b>pH stability</b>	3-13	3-13
<b>Temperature</b>	4-30°C	4-30°C
<b>Chemical stability</b>	All commonly used aqueous buffers, 1 M sodium hydroxide (NaOH), 6 M guanidine hydrochloride, 30% isopropanol, and 70% ethanol.	
<b>Storage</b>	20% ethanol at 4°C to 25 °C	



# Instruments & Parts

GALAK provides selected instruments and parts that are used in HPLC systems.

- Packing system for HPLC columns
- High-pressure Injection Pump
- Low-pressure Glass Column
- Injection Loop
- Oligo Synthesis Column
- Empty HPLC Column
- Accessories for HPLC system



# Packing System For HPLC Column

GLK Packing Systems are designed for packing analysis, semi-preparative and preparative HPLC columns.

GLK 1000, designed for packing analytical columns only, is suitable for the packing of conventional silica-gel and polymer HPLC columns.

GLK 2000, with higher pressure and power, are designed for both analytical and preparative columns with inner diameter 10~50mm.

Customized homogenate tanks suitable for HPLC columns are provided optional.

GLK HPLC Column Packing System is widely used in many famous universities and research institutions such like Tsinghua University, Sichuan University, Zhengzhou University, Dalian Institute of Chemical Physics Dalian Ocean University.

## Advantage:

- Desktop level device. It could be put in a moving cart.
- Compressed air for power. No electricity required.
- The main machine is practical, robust and durable.
- Customized accessories such as homogenate tanks and connectors.



## Service:

- One year warranty
- Free replacement parts
- Free online training for operation and maintenance
- Recovery of old equipment

## Parameters:

	GLK1000	GLK2000
Column ID	2.0/2.1/4.6/10 mm	4.6/10/20/30/50 mm
Output Pressure	9800 psi	15000 psi
Flow Rate	3.3L/min	3.3L/min
Output Power	1.5hp	2hp
Air Cylinder	Single	Double
Size	50*40*20cm	60*40*20cm
Weight	15kg	20kg

## Control Panel

- 1 Pressure gauge
- 2 Pressure regulator
- 3 Liquid inlet:
- 4 Inlet A:
- 5 Inlet B:
- 6 Liquid outlets:



## Homogenate Tank & Connector



## Equipment Accessories

Standard Parts	Optional Parts
Operation instruction	Air compressor
Pneumatic booster pump	Air purification system
Control panel	Homogenate tanks
Homogenate tank support (mini)	Column connection (ID 2.1-50mm)
Stainless steel connections	Empty HPLC column (ID 2.1-50mm)
Stainless steel column	Packing materials

Notice:

- 1. Nitrogen compressed air cylinder could replace air compressor and air purification system.
- 2. The compressed air must be purified before entering the packing system.

# High-pressure Injection Pump

## Eldex Optos Electronic Metering Pump

Eldex pumps are produced in the USA. Integrating the latest technology and electronics, Eldex's Optos Series designed and manufactured with reciprocating piston pumps for a wide variety of applications.

- Smooth fluid delivery from intra
- RPM stepper motor control
- Remote flow rate control via analog signal or RS232 command
- Optional liquid ends for expanded flow rate range
- Integrated piston wash system

### With upgrade to Plus Version

- Pressure monitoring with high and low pressure limits
- Integrated low volume pulse damper



## Model 1

	Flow Rate (mL/min)	Max. Pressure (psi)	Piston Diameter (in.)	Piston Stroke (in.)	Model
<b>316 stainless steel</b>	0.002 - 2.5	6000	3/32	.125	1LM
	0.003 - 5	6000	1/8	.125	1SM
	0.01 - 20	3000	1/4	.125	1HM
	Flow Rate (mL/min)	Max. Pressure (psi)	Piston Diameter (in.)	Piston Stroke (in.)	Model
<b>PEEK</b>	0.002 - 2.5	4000	3/32	.125	1LI
	0.003 - 5	4000	1/8	.125	1SI
	0.01 - 20	3000	1/4	.125	1HI

### Model 2

316 stainless steel	Flow Rate (mL/min)	Max. Pressure (psi)	Piston Diameter (in.)	Piston Stroke (in.)	Model
	0.003 - 5	6000	3/32	.250	2LM
	0.01 - 10	6000	1/8	.250	2SM
	0.02 - 40	1500	1/4	.250	2HM
PEEK	Flow Rate (mL/min)	Max. Pressure (psi)	Piston Diameter (in.)	Piston Stroke (in.)	Model
	0.003 - 5	4000	3/32	.250	2LI
	0.01 - 10	4000	1/8	.250	2SI
	0.02 - 40	1500	1/4	.250	2HI

### Model 3

316 stainless steel	Flow Rate (mL/min)	Max. Pressure (psi)	Piston Diameter (in.)	Piston Stroke (in.)	Model
	0.01 - 10	3000	3/32	.500	3LM
	0.01 - 20	1500	1/8	.500	3SM
	0.04 - 80	750	1/4	.500	3HM
PEEK	Flow Rate (mL/min)	Max. Pressure (psi)	Piston Diameter (in.)	Piston Stroke (in.)	Model
	0.01 - 10	3000	3/32	.500	3LI
	0.01 - 20	1500	1/8	.500	3SI
	0.04 - 80	750	1/4	.500	3HI

### Optos Plus Model: Minimize Pulsation, Monitor Pressure

Add Plus to your Optos Series pump to integrate a pulse damper to further reduce pulsation and have the ability to monitor pressure and set high and low pressure limits. Plus is available on L and S piston pumps.

	Flow Rate* (mL/min)	Max. Pressure (psi)	Piston Diame- ter (in.)	Piston Stroke (in.)	Model
316 stainless steel	0.002 - 2.5	6000	3/32	.125	1LMP
	0.003 - 5	6000	1/8	.125	1SMP
	Flow Rate* (mL/min)	Max. Pressure (psi)	Piston Diame- ter (in.)	Piston Stroke (in.)	Model
PEEK	0.002 - 2.5	4000	3/32	.125	1LIP
	0.003 - 5	4000	1/8	.125	1SIP

### Common Specifications

**Wetted Parts:** For 316: Type 316 stainless steel, sapphire, ruby, gold, UHMW Polyethylene, CTFE (consult factory for other options)

For PEEK: PEEK, sapphire, ruby, UHMW Polyethylene, inert polymers

**Reproducibility:** typically +/-0.3%

**Viscosity Limit:** 500 centipoise

**Tubing Connections:** Consult factory for additional options

For 316: Low and Standard Flow Liquid End pumps:

Inlet valves: 1/4"-28 plastic fitting for 1/8" plastic tubing Outlet valves: 10-32 tube nut, ferrule for 1/16" tubing

For 316: High Flow Liquid End pumps: Inlet and outlet valves: 1/8" Swagelok®

For PEEK: Low and Standard Flow Liquid End pumps:

Inlet valves: 1/4"-28 plastic fitting for 1/8" plastic tubing Outlet valves: 1/4"-28 plastic fitting for 1/16" tubing

For PEEK: High Flow Liquid End pumps: Inlet and outlet valves: 1/4"-28 plastic fitting for 1/8" tubing

**Control Options:** Voltage: 0-5 VDC; Current Loop: 4-20 mA; RS232

**Contact Closure & Outputs:** Remote run and stop, pressure output (0-5V), pump error output.

**Dimensions:** 10 cm W x 23 cm H x 24 cm D (4" W x 9" H x 9.5" D)

**Weight:** with Plus: 6.4 kg (14 lbs.); w/o Plus: 5.3 kg (11.75 lbs.)

**Electrical:** 100-230VAC (+/- 10%); 50/60Hz

**VA Rating:** 80

**CE Certified**

# Single-layer Glass Columns


GALAK chromatography glass columns are designed for the standard liquid chromatography of macromolecules. The design is based on high reproducibility and precision results. The columns are compatible with aqueous solution and organic solvent in liquid chromatography. The GALAK column head is applied with chromatography column which includes a head and a bottom. A larger scale of column bed height can be obtained if the bottom part is replaced with another GALAK column head.

## Advantages

- Pressure-resistant borosilicate glass, visualization and stability
- Supporting foot, adjustable level, convenient for users to use
- Reasonable price, high cost performance
- Reproducibility, excellent column efficiency and reliable results
- Zero dead volume structural connections



No	Name	Material
1	Glass Tube	Borosilicate glass
2	Screw	POM
3	Compression bar	POM
4	Screw cap	POM
5	Piston	PP
6	Compression nut	POM
7	Screw cap	POM
8	Filter	PA/PP
9	Support net	PP
10	Screw cap	POM
11	O-ring	VITON (EPDM)
12	Compression bar	POM
13	Locking ring	ETFE
14	Connector	PEEK
15	Capillary	PFA
16	Support Legs	304





<b>Working Temperature</b>	4-40°C
<b>pH Range</b>	1-14
<b>Chemical Stability</b>	Tolerant to salt, acid, alkali, and a small number of organic solvents alcohols, ketones, phenols.
<b>Column Material</b>	Borosilicate glass
<b>Column Head Material</b>	PTFE
<b>Thread-end Material</b>	PEEK
<b>Seal Ring Material</b>	PTFE/EPDM
<b>Tubing Material</b>	1/16&1/8
<b>Connector Material</b>	PEEK 1/16&1/8

No.	Internal Diameter (mm)	Length (mm)	One-side Adjustable Type		Double-side Adjustable Type		Pressure (bar)
			Volume (mL)	Bed Height (cm)	Volume (mL)	Bed Height (cm)	
YS16/200	16	200	4-30	2-14.5	0-30	0-14.5	7
YS16/400	16	400	46-72	22-34.5	17-72	8.5-34.5	7
YS16/1000	16	1000	173-199	82-94.5	144-199	68.5-94.5	7
YS26/200	26	200	10-73	2-14.5	0-73	0-14.5	7
YS26/400	26	400	111-174	22-34.5	43-174	8.5-34.5	7
YS26/1000	26	1000	415-479	82-94.5	347-479	68.5-94.5	7

# BSXK Double-layer Glass Columns

BSXK glass columns are made of borosilicate glass. They allow visual inspection of media bed and exhibit excellent chemical resistance. Column packing can be performed using either a packing reservoir or extra column tube attached with a packing connector. QuickLock of the adapter shaft facilitates rapid and easy movement of the adapter, simplifying adjustments of the bed height and cleaning. Adapter plunger gives a uniform flow which maintains the integrity of the packed bed during operation.

## Advantages

- Pressure resistant glass column tube can be used to clearly observe the column bed
- The column head ensures uniform buffer distribution and repeated loading
- Unique support foot design can adjust the level after placement
- Various specifications, inner diameter: 10, 16, 26, 50mm
- Maximum compression at one end is 17cm, and 2 \* 17cm at both ends
- Direct connection with mainstream purification system



Number	Name	Number	Name
1	Column Head	14	Support net
2	Column tube	15	Sieve
3	Column base	16	Locking pressure ring
4	Pipe joint	17	Gasket
5	Capillary	18	Sealing ring
6	Compression nut	19	Connector
7	Screw cap	20	O-ring
8	Screw	21	Acrylic tube
9	Pressure bar	22	Glass tube
10	Seal lock ring	23	Lower piston
11	Piston gland	24	Lower screw cap
12	Upper piston	25	Support leg
13	O-ring		



<b>Working Temperature</b>	4-40°C
<b>pH Range</b>	1-14
<b>Chemical Stability</b>	Tolerant to salt, acid, alkali, and a small number of organic solvents alcohols, ketones, phenols.
<b>Column Material</b>	Borosilicate glass
<b>Column Head Material</b>	PTFE
<b>Thread-end Material</b>	PEEK
<b>Seal Ring Material</b>	PTFE/EPDM
<b>Tubing Material</b>	1/16&1/8
<b>Connector Material</b>	PEEK 1/16&1/8
<b>Max. Pressure</b>	5 bar

No.	Internal Diameter (mm)	Length (mm)	One-side Adjustable Type		Double-side Adjustable Type	
			Volume (mL)	Bed Height (cm)	Volume (mL)	Bed Height (cm)
BSXK10/150	10	150	-		0-12	0-11
BSXK10/200	10	200	-		0-25	0-16
BSXK10/400	10	400	-		0-68	21-36
BSXK10/1000	10	1000	-		0-150	81-96
BSXK16/200	16	200	4-30	2-14.5	0-30	0-14.5
BSXK16/400	16	400	46-72	22-34.5	17-72	8.5-34.5
BSXK16/700	16	700	109-136	52-64.5	81-136	38.5-64.5
BSXK16/1000	16	1000	173-199	82-94.5	144-199	68.5-94.5
BSXK26/200	26	200	10-73	2-14.5	0-73	0-14.5
BSXK26/400	26	400	111-174	22-34.5	43-174	8.5-34.5
BSXK26/700	26	700	263-326	54-64.5	195-326	38.5-64.5
BSXK26/1000	26	1000	415-479	82-94.5	347-479	68.5-94.5
BSXK50/200	50	200	19-275	1-14	0-275	0-14
BSXK50/300	50	300	215-471	11-24	0-471	0-24
BSXK50/600	50	600	804-1060	41-54	549-1060	28-54
BSXK50/1000	50	1000	1589-1849	81-94	1334-845	68-94

# Low-Pressure Chromatography Column

Low-pressure chromatography columns are pressure compressible glass columns designed for hygienic operation and simple, efficient loading, primarily for process development or biopharmaceutical production.

## Patent column head sealing technology

- The lever-pressing sealing structure was used with high reliability, which prevents the problem that the pneumatic mechanism easy to leak and invalid.
- The expansion structure of pressurizing-down style gasket ring prevents column head departing from bed caused by the pull-up structure.
- Minimized Hold-up Volumes, Easy to clean and change the seal.



## Patent column head rotating structure

- The column head rotates by the rotating screw of the column pipe, which is on the upper surface of the flange plate. After rotating in place, the second screw needs to be inserted. Media packing can be done after rotating the column head. It is easy to operate, without carrying out the column head.

## Predictable linear scale-up

- Fix condition: Linear flow rate, buffer, packing material, bed height, sample concentration, and pH, sample volume, and bed volume ratio.
- Scale-up condition: Column I.D., volume flow rate, sample volume.

Product	Column Inner Diameter (mm)	Sectional Area (cm <sup>2</sup> )	Column Height (mm)	Column Bed Height (cm)		Column Bed Volume (L)		Max. Pressure (bar)	Net Weight (Kg)
				Min	Max	Min	Max		
MPC100/500	70	38.5	500	0	35	0	1.4	8	14
MPC100/750	70	38.5	950	40	80	1.5	3.1	8	14
MPC100/500	100	78.5	500	0	35	0	2.7	8	18
MPC100/750	100	78.5	750	20	60	1.6	4.7	8	20
MPC100/950	100	78.5	950	40	80	3.1	6.3	8	21
MPC140/500	140	154	500	0	35	0	5.4	6	30
MPC140/750	140	154	750	20	60	3.1	9.2	6	33
MPC140/950	140	154	950	40	80	6.2	12.3	6	35
MPC200/500	200	314	500	0	35	0	11	6	36
MPC200/750	200	314	750	20	60	6.3	18.8	6	39
MPC200/950	200	314	950	40	80	12.6	25.1	6	42
MPC300/500	300	706.5	500	0	35	0	24.7	4	58
MPC300/750	300	706.5	750	20	60	14.1	42.4	4	63
MPC300/950	300	706.5	95	40	80	28.2	56.5	4	67
MPC450/500	450	1560	50	0	35	0	55.6	3	230

## Adaptor For Glass Column

The adaptor is matched with the glass column.

YS columns: 130ml / 320ml

BSXK columns: 30ml / 130ml / 320ml / 420ml / 1200ml

HT columns: HT10/110, HT26/100, HT49/100



## Injection Loop

GALAK injection loop can be incorporated into a pressurized packing device for large volume samples and used in conjunction with the sampling valve.

### Advantage

- Flexible injection of different volumes of samples
- High reproducibility and recovery
- Internal dynamic seal allows the sample to be released

### Sample Volume

- 10mL
- 50mL
- 150mL



# Oligo Synthesis Column

Adjustable Oligo synthesis columns are designed for oligonucleotide synthesis and can withstand the harsh organic conditions of oligonucleotide synthesis.

- Adjustable column bed height recommended working height between 3 and 10 cm
- Column diameter of 35 mm and volume range of 10 mL to 100 mL
- The column can withstand the harsh organic conditions in synthesis
- The columns are easy to handle and the solid phase carriers are easily packaged

Oligo columns are also available for large-scale synthesis in 70mm, 100mm, 200mm and 350mm diameters.



Product	Inner Diameter	Height
Oligo35 Column	35mm	150mm
Oligo50 Column	50mm	250mm
Oligo100 Column	100mm	500mm

Small stainless steel synthesis columns are designed as fixed volume synthesis column reactors (equipped with filters and seals) for oligonucleotide synthesis.

- Synthetic columns are manufactured to high standards to withstand the harsh organic conditions of oligonucleotide synthesis
- Made of 316L stainless steel
- Available in 1.2 ml, 6.3 ml, 12 ml, 24 ml and 48 ml sizes

Oligo columns are also available for large-scale synthesis in 70mm, 100mm, 200mm and 350mm diameters.



Product	Volume	Inner Diameter	Height
Oligo1.2	1.2ml	10mm	15mm
Oligo6.3	6.3ml	20mm	20mm
Oligo12	12ml	27mm	21mm
Oligo24	24ml	35mm	25mm
Oligo48	48ml	44mm	32mm



## HPLC Column Hardware

GALAK empty columns are designed for HPLC system. The column tube is made of high quality 316L stainless steel tube and precision machining.

- high-purity 316 L stainless steel tubing material.
- Passivated column tube for good acid and alkali resistance..
- Roughness of tube inner wall Ra0.3.
- Circular cross-section.
- High consistency and reproducibility.
- Sieve plates are easy to replace and clean.
- Low dead volume.
- OEM is available.



### Column Size:

- Inner diameter: 2.1mm, 3.0mm, 4.0mm, 4.6mm, 7.8mm, 10mm, 20mm, 21.2mm, 30mm, 50mm
- Length: 25mm, 30mm, 50mm, 100mm, 150mm, 250mm, 300mm, 500mm
- Frits: 2um (optional 0.5um, 1um, 3um, 5um, 10um)
- Connector: 1/16" (optional 1/8")



## Frits For HPLC column

Inner diameter: 2.1mm, 3.0mm, 4.0mm, 4.6mm

Material: 316L stainless steel & PEEK



## In-filter for HPLC column

### Type:

- 10mm I.D.
- 20mm I.D.
- 30mm I.D.



# **GALAK**

## CHROMATOGRAPHY

---

GALAK Chromatography Technology Co., LTD.

88th Meiliang Rd, Mashan, Wuxi, China

Tel: +86 510 85992929 Fax: +86 510 85104949

Email: [sales@galak-tech.com](mailto:sales@galak-tech.com)

[www.galak-tech.com\(CN\)](http://www.galak-tech.com(CN)) [www.galaklc.com\(US\)](http://www.galaklc.com(US))